

Laboratory Sample Condenser



Headers details

This LSC unit ranges have been designed to allow clean steam and Water For Injection (WFI) samples to be taken quickly and easily whilst providing a sterile environment for testing and can be operated with mains water as the cooling medium. Availability of aseptic sample valve allow fine control of sample flow during testing. The capacity of the LSC range is dependent on the system pressure and temperature, together with the temperature and flow rate of the cooling water supply.

		Aerre Inox s.r.l.	
L'Iniziazione Cesare - 20090 FIESCO (CR) Italy Tel. 0374-370038 - Fax 0374-370633			
TYPE - TIPO		SERVICE - SERVIZIO	
ITEM	SHOP No. / No. FABR.	YEAR - ANNO	
	SHELL - MANTELLO	TUBES - TUBI	
CAPACITY CAPACITÀ			
DESIGN PRESSURE PRESSIONE DI PROGETTO	BAR		
DESIGN TEMPERATURE TEMPERATURA DI PROGETTO	°C		
HYDR. TEST PRESSURE PRESS. PROVA IDRAULICA	BAR		
FLUID GROUP GRUPPO FLUIDI			

Clamped Headers

